

## **SUPPLEMENTARY MATERIAL**

If not otherwise stated, all chemicals were from Sigma-Aldrich (Merck, Darmstadt, Germany).

### **Preparation of fecal extracts**

Details can be found in the supplementary material. Frozen fecal samples were grinded into fine powder by cryomilling in liquid nitrogen using a fecal-specific milling program (SamplePrep 6875D, SPEX, Metuchen, USA) [1]. Pulverized samples were used for the preparation of protein extracts but also the extraction of DNA. Powdered fecal samples were resuspended in lysis buffer (50 mM Tris-HCl, 150 mM NaCl, pH 8), 300 mg sample in 1.5 mL buffer and were subjected to bead-beating at 30 Hz for 2 x 10 min (Stainless Steel Beads 5 mm, Qiagen, Hilden, Germany); Mixer Mill MM400, Retsch, Haan, Germany). After overnight rotation at 4°C, samples were centrifuged at 14,770 x g for 45 min at 4°C. The resulting supernatant was used as a first fraction. The pellet was further extracted using lysis buffer with 4 M urea for protein stabilization and repeating the extraction steps [2]. This second fraction was combined with the first and used for analysis. A Bradford assay was used to determine the protein concentration of each fecal protein extract. 2 µL extract was diluted in 798 µL H<sub>2</sub>O and 200 µL assay dye was added (Protein assay dye 5000006, BIO-RAD, Winninglaan, Temse, Belgium). The quantification was performed according to the manufacturers' protocol, using a standard curve of known concentrations of bovine serum albumin (BSA), absorption was measured at 595 nm.

### **Enzyme-linked immunosorbent assay (ELISA)**

Also, fecal Calprotectin (fCal) was assessed as a measure for intestinal inflammation. Ground frozen fecal samples of peanut-allergic children and healthy age-matched controls were used according to the manufacturer's instructions provided in the kit (EK-CAL, Bühlmann Laboratories AG, Schönenbuch, Switzerland). The established cut-off for fCal is 50 µg/g, with levels higher than that indicating intestinal inflammation. In our cohort, fCal results did not discriminate between peanut-allergic children and controls. Healthy children had a mean fCal value of 75.8 µg/g (range 11.4-260.9 µg/g of fecal pellet), while peanut-allergic patients had mean value of 32.6 µg/g (range 0-185.6 µg/g). Indeed, it has been shown before that healthy children may have fCal levels over this cut-off level and that levels decrease with age [3,4].

### **Multiplex analysis for deep IgE-profiling**

Reactivity profiles of specific IgE (sIgE) antibodies in participant sera and fecal extracts were assessed using the ALEX2 macroarray with 298 allergens (MADX, Vienna, Austria) [5,6]. The assay was performed according to the manufacturer's protocol. Briefly, 100  $\mu$ L of 1:5-diluted fecal extract or serum samples were used. The readout was done using array-specific software. Per patient, fecal sIgE values above the cut-off for positivity (0.3 kU<sub>A</sub>/L) were added up to yield the sum of sIgE (e.g. Table 1). The mean sum of fecal sIgE was 31.2 kU/L (range 5.22-75.82 kU/L).

**REFERENCES:**

1. Hickl O, Heintz-Buschart A, Trautwein-Schult A, Hercog R, Bork P, Wilmes P, et al. Sample Preservation and Storage Significantly Impact Taxonomic and Functional Profiles in Metaproteomics Studies of the Human Gut Microbiome. *Microorganisms*. 2019;7(9):367.
2. Kuehn A, Scheuermann T, Hilger C, Hentges F. Important variations in parvalbumin content in common fish species: a factor possibly contributing to variable allergenicity. *Int Arch Allergy Immunol*. 2010;153(4):359-66.
3. Velasco Rodríguez-Belvís M, Viada Bris JF, Plata Fernández C, García-Salido A, Asensio Antón J, Domínguez Ortega G, et al. Normal fecal calprotectin levels in healthy children are higher than in adults and decrease with age. *Paediatr Child Health*. 2020;25(5):286-292.
4. Roca M, Rodríguez Varela A, Carvajal E, Donat E, Cano F, Armisen A, et al. Fecal calprotectin in healthy children aged 4-16 years. *Sci Rep*. 2020;10(1):20565.
5. Heffler E, Puggioni F, Peveri S, Montagni M, Canonica GW, Melioli G. Extended IgE profile based on an allergen macroarray: a novel tool for precision medicine in allergy diagnosis. *World Allergy Organ J*. 2018;11(1):7.
6. Czolk R, Ruiz- Castell M, Hunewald O, Wanniang N, Le Coroller G, Hilger C, et al. Novel, computational IgE- clustering in a population- based cross- sectional study: Mapping the allergy burden. *Clin Transl Allergy*. 2023;13(9):e12292.